

Wir begleiten  
Ihre erfolgreiche  
Getränkeherstellung

**SCHLIESSMANN  
SCHWÄBISCH HALL**



Tel. 07 91 - 9 71 91-0 • Fax 9 71 91-25  
C. Schliessmann Kellerei-Chemie GmbH & Co.KG  
Auwiesenstr. 5 • D-74523 Schwäbisch Hall

Getränkeanalytik

# BENTOTEST<sup>®</sup> according to Dr. L. Jakob

version 04/2007

- rapid determination of the amount of bentonite needed  
for wine and juices -

page 1/2

- BENTOTEST<sup>®</sup>-solution (Nr. 2508) for colorintensive red wines for red wine
- BENTOTEST<sup>®</sup>- solution (Nr. 2509) or strongly neutralized white wines up to bright red wines for white wine (colorless)
- BENTOTEST<sup>®</sup>- solution (Nr. 2510) for white wine for white wine (yellow)

## Preliminary test:

To find out whether bentonite fining is necessary for the drink under investigation, the following preliminary test is carried out:

1 part of BENTOTEST<sup>®</sup> reagent is added to 10 parts of filtered wine (room temperature!). The addition of 1 part of reagent to 10 parts of drink need only be approximate. In practice, it is easiest to start from the total amount of drink in the flask (about 50cm<sup>3</sup>) and to add about 5cm<sup>3</sup> of the BENTOTEST<sup>®</sup> reagent by means of the measuring beaker provided. In wines in need of bentonite treatment, the drink turns turbid. With a little practice, the amount of bentonite needed can be concluded from the degree of turbidity.

## The rule of thumb is as follows:

Slight turbidity:	50 to 100g of bentonite per hectolitre of wine
Medium turbidity:	100 to 250g of bentonite per hectolitre of wine
Extreme turbidity:	250 to 400g of bentonite per hectolitre of wine

To determine the exact amount of bentonite, fining experiments have to be carried out as follows.

## Preliminary fining experiment:


1. Fill the drink from the cellars into the flask up to the O mark.
2. Adjust exactly to the O mark by means of the pipette provided.
3. Vigorously shake the bottle with the bentonite suspension (yellow label).
4. Add the well-mixed bentonite suspension to the drink in the flask in the amount estimated in the preliminary test.
5. Close the flask by placing a thumb over it, and shake for 2-3 minutes.
6. Then filter through BENTOTEST<sup>®</sup> special folded filters in an Erlenmeyer flask.
7. Add 1 part of BENTOTEST<sup>®</sup> reagent to 10 parts of filtered wine (follow the instructions in the preliminary test).

**The wine is protein-free when no turbidity appears.**

The amount of bentonite needed is obtained most accurately by running 3 fining experiments for each wine or fruit juice with rising amounts of bentonite, and observing the amount which just makes the drink protein-free.

It should be stressed that the bentonite suspension used should always be from the bentonite employed in practical fining in the cellars.

## BENTOTEST<sup>®</sup> equipment and replacement parts



**A quick method to determine the Bentonite requirements needed in the fining of wines and juices.**

**It gives guarantee**

through absolutely reliable evidence of the presence of protein, whereby all protein species are comprehended.

**It helps in the economizing Bentonite**

through accurate determination of the required quantity of product needed for the fining.

**Fast review**

ensuing the fining and before the filling, which is of greatest importance, especially for warm filling.

BENTOTEST<sup>®</sup> equipment complete in plastic case

(Nr. 2500)

**includes:**

2 Bentotest flasks, 3 Erlenmeyer flasks, 3 plastic hoppers, 1 plastic measuring beaker, 1 pipette, 50 special paper filters, 1 bottle Bentotest-solution, 1 bottle bentonite-suspension, directions for use, all include in a plastic case

**Replacement parts:**

Bentotest flask	(Nr. 2502)
Erlenmeyer flask	(Nr. 2503)
Plastic hopper	(Nr. 9950)
Plastic measuring beaker	(Nr. 2504)
Pipette	(Nr. 2505)
pack special paper filter (100 piece)	(Nr. 2506)
BENTOTEST <sup>®</sup> -solution for white wine (colorless), 250 ml	(Nr. 2509)
BENTOTEST <sup>®</sup> -solution for white wine (yellow), 250 ml	(Nr. 2510)
BENTOTEST <sup>®</sup> -solution for red wine, 250 ml	(Nr. 2508)
Ca-Bentonite-suspension 250 ml	(Nr. 2511)
Na-Ca-Bentonite-suspension 250 ml	(Nr. 2512)